

POWER MOBILITY WHEELCHAIR

ESCAPE HP1

A steel-framed model serves as an economical alternative.



ESCAPE DX HP5

An Aluminum-framed model features Aluminum die-casting lugs instead of welded ones.



ESCAPE LX HP8



Adjustable Backrest
(90~110 degrees)

- HP8 Escape LX is a strong, durable electric wheelchair ideal for portability and general indoor/outdoor use
- Available in 16", 18", and 20" seat widths. Adjustable backrest



ESCAPE S P12S

A model incorporates a rear-wheel suspension system, and thus provides great comfort and ensures a good ride.



ESCAPE SX P12SX

A strong and enduring model, with 9" front castors and 14" rear wheels installed, is ready for a faster stroll.



ESCAPE AX P21

- Detachable and height adjustable armrest
- Height adjustable backrest at most 10cm
- Firm backrest support bar



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
HP1	135kgs (300lbs)	1090mm x 630mm	62kgs (136lbs)	200W x 2pcs	6KPH (3.8MPH)	12V 36Ah x 2pcs	30km (20miles)	330mm x 50mm (13" x 2")	5 Amp Off Board	865mm (34")
HP5	135kgs (300lbs)	1100mm x 685mm	62kgs (136lbs)	200W x 2pcs	6KPH (3.8MPH)	12V 36Ah x 2pcs	30km (20miles)	320mm x 60mm (12.5" x 2.375")	5 Amp Off Board	865mm (34")
HP8	135kgs (300lbs)	1085mm x 640mm	66kgs (145lbs)	200W x 2pcs	6KPH (3.8MPH)	12V 36Ah x 2pcs	30km (20miles)	320mm x 60mm (12.5" x 2.375")	5 Amp Off Board	865mm (34")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P12S	135kgs (300lbs)	1120mm x 610mm	73kgs (161lbs)	200W x 2pcs	6KPH (3.8MPH)	12V 36Ah/50Ah x 2pcs	30km (20miles) / 45km (28miles)	320mm x 60mm (12.5" x 2.375")	5 Amp Off Board	900mm (35.5")
P12SX	160kgs (350lbs)	1130mm x 660mm	83kgs (183lbs)	600W x 2pcs	12KPH (7.5MPH)	12V 36Ah/50Ah x 2pcs	30km (20miles) / 45km (28miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	950mm (37")
P21	135kgs (300lbs)	1030mm x 630mm	80kgs (176lbs)	450W x 2pcs	6KPH (3.8MPH)	12V 50Ah x 2pcs	45km (28miles)	320mm x 60mm (12.5" x 2.375")	5 Amp Off Board	850mm (33.5")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

POWER MOBILITY WHEELCHAIR

ESCAPE HD HP1HD




A heavy-duty steel-framed wheelchair, structured with a double-cross brace design and equipped with 4-pole high power motors, offers seat width options of 22" and 24" for the heavyweights.



ESCAPE LXHD HP8HD



HP8HD is a new addition to Heartway's renown Bariatric range of products. Its double-cross bar design greatly strengthens the structure of the product. It offers seat width options of 24", 26" and 28". Also as a folding chair, at the same time it offers high portability for our big friends who don't mind travelling a little bit!



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
HP1HD	180kgs (400lbs)	1150mm x 750mm	82kgs (180lbs)	700W 4-Pole Motor x 2pcs	10KPH (6.25MPH)	12V 50Ah x 2pcs	20km (13miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	1000mm (39")
HP8HD	225kgs (500lbs)	1170mm x 910mm	85kgs (187lbs)	700W 4-Pole Motor x 2pcs	10KPH (6.25MPH)	12V 62Ah x 2pcs	30km (20miles)	360mm x 75mm (14" x 3")	8 Amp Off Board	1140mm (45")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

MAXX P3D



A mid-wheel drive power wheelchair, coming with a unique independent suspension system, offers great stability and lots of comfort and is featured with an obstacle height of 6cm and more.



MAXX R P3DR



A mid-wheel drive power wheelchair with a rehab seating system provides a wonderful driving experience due to its very special independent suspension system.



MAXX RT P3DRT



A mid-wheel drive model featuring a tilt-in-space rehab seating system offers great mobility not just with the chair but also on the chair.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P3D	160kgs (350lbs)	1080mm x 680mm	107kgs (235lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	30km (20miles)	305mm x 100mm (12" x 4")	5 Amp Off Board	550mm (21.625")
P3DR	160kgs (350lbs)	1080mm x 680mm	107kgs (235lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	30km (20miles)	330mm x 100mm (13" x 4")	5 Amp Off Board	550mm (21.625")
P3DRT	160kgs (350lbs)	1080mm x 680mm	107kgs (235lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	30km (20miles)	330mm x 100mm (13" x 4")	5 Amp Off Board	550mm (21.625")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

POWER MOBILITY WHEELCHAIR



RUMBA HP3

A strong and enduring model, coming with 9.5" front castors and 14" rear wheels, is great for a faster stroll and offers sufficient obstacle-climbing height.



RUMBA S HP4

A light and agile model that is just perfect for home care.



RUMBA SF P4F

A front-wheel drive model which has great agility and a small turning radius.



RUMBA SR P4R

A rear-wheel drive model that has easy maneuverability and offers a good ride.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
HP3	135kgs (300lbs)	1000mm x 610mm	98kgs (216lbs)	450W x 2pcs	8KPH (5MPH)	12V 36Ah/50Ah x 2pcs	30km (20miles) / 45km (28miles)	330mm x 100mm (13" x 4")	5 Amp Off Board	825mm (32.5")
HP4	135kgs (300lbs)	1100mm x 630mm	83kgs (183lbs)	450W x 2pcs	7KPH (4.4MPH)	12V 36Ah/50Ah x 2pcs	30km (20miles) / 45km (28miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	785mm (31")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P4F	135kgs (300lbs)	1110mm x 610mm	93kgs (205lbs)	450W x 2pcs	7KPH (4.4MPH)	12V 36Ah/50Ah x 2pcs	30km (20miles) / 45km (28miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	530mm (20.875")
P4R	135kgs (300lbs)	1040mm x 610mm	93kgs(205lbs)	450W x 2pcs	7KPH (4.4MPH)	12V 36Ah/50Ah x 2pcs	30km (20miles) / 45km (28miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	755mm(29.5")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

POWER MOBILITY WHEELCHAIR



TIARA P4AS

A model gives a clean look and highlights a simple but functional design.



TANGO TC1

A model with an attendant-control mechanism makes it effortless to transit patients. Breathable and patent-protected material is used on the seating.




TIARA R P4ASR

Tiara with a rehab seating system provides more width and depth adjustments.



TANGO S TC1S

- Equipped with nylon back-rest cushion and nylon seat cushion
- Detachable armrest with height adjustment available
- Attendant wheelchair – TC1S provides the most support and comfort for users



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P4AS	135kgs (300lbs)	1100mm x 660mm	90kgs (198lbs)	450W x 2pcs	7KPH (4.4MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	790mm (31")
P4ASR	135kgs (300lbs)	1100mm x 660mm	94kgs (207lbs)	450W x 2pcs	7KPH (4.4MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	860mm (34")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
TC1	135kgs (300lbs)	1160mm x 590mm	77kgs (169lbs)	450W x 2pcs	6KPH (3.8MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	750mm (29.5")
TC1S	135kgs (300lbs)	1160mm x 590mm	77kgs (169lbs)	450W x 2pcs	6KPH (3.8MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	750mm (29.5")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

POWER MOBILITY WHEELCHAIR

ALLURE

HP6

A mid-wheel drive model offers a small turning radius and presents a classic look.



ALLURE R

HP6R

A rehab seating system incorporated provides many seat width and seat depth options.



ALLURE S

HP6S

Equipped with suspension system.



ALLURE RT

HP6RT

A tilt-in-space rehab seating system provides the comfort you've never experienced before!



SAHARA R

HP7R

A model featuring front, middle, and rear independent suspension systems, equipped with a rehab seating system highlighting the design of Ergonomics technology.



SAHARA RT

HP7RT

Incorporated with a tilt-in-space rehab seating system and featured with quality linear guide ways, is a great example of applying Ergonomics technology in power mobility products.



SAHARA K

HP7K

A model features a comfortable captain seat.



SAHARA KX

HP7KX

Available with 6.5" wider rear wheels.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
HP6	160kgs (350lbs)	1110mm x 620mm	97kgs (213lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	40km (25miles)	330mm x 100mm (13" x 4")	5 Amp Off Board	660mm (26")
HP6S	160kgs (350lbs)	1110mm x 620mm	97kgs (213lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	40km (25miles)	330mm x 100mm (13" x 4")	5 Amp Off Board	660mm (26")
HP6R	160kgs (350lbs)	1110mm x 620mm	101kgs (222lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	40km (25miles)	330mm x 100mm (13" x 4")	5 Amp Off Board	660mm (26")
HP6RT	160kgs (350lbs)	1110mm x 620mm	114kgs (251lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	40km (25miles)	330mm x 100mm (13" x 4")	5 Amp Off Board	660mm (26")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
HP7R	205kgs (450lbs)	1170mm x 635mm	113kgs (249lbs)	700W 4-Pole Motor x 2pcs	10KPH (6.25MPH)	12V 50Ah/62Ah x 2pcs	20km (13miles) / 30km (20miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	865mm (34")
HP7RT	205kgs (450lbs)	1170mm x 635mm	135kgs (300lbs)	700W 4-Pole Motor x 2pcs	10KPH (6.25MPH)	12V 50Ah/62Ah x 2pcs	20km (13miles) / 30km (20miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	865mm (34")
HP7K	205kgs (450lbs)	1170mm x 830mm	128kgs (282lbs)	700W 4-Pole Motor x 2pcs	10KPH (6.25MPH)	12V 50Ah/62Ah x 2pcs	20km (13miles) / 30km (20miles)	360mm x 160mm (14" x 6.5")	5 Amp Off Board	865mm (34")
HP7KX	205kgs (450lbs)	1170mm x 830mm	128kgs (282lbs)	700W 4-Pole Motor x 2pcs	10KPH (6.25MPH)	12V 50Ah/62Ah x 2pcs	20km (13miles) / 30km (20miles)	360mm x 160mm (14" x 6.5")	5 Amp Off Board	865mm (34")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

POWER MOBILITY WHEELCHAIR

VISION C P13C

A model with a captain seat incorporated gives an elegant look.



VISION P13

A model which highlights an elevating power actuator and a double cross brace design is indeed classic and elegant. The seat, with great stability and quietness, elevates by 30 cm, which makes it rather easy for the users to grab articles at higher places.



VISION R P13R

A model with a rehab seating system incorporated gives an elegant look.



VISION RT P13RT



A model comes with a tilt-in-space rehab seating system and features linear guideways of high stiffness, precise accuracy, and great loading capacity.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P13C	160kgs (350lbs)	1120mm x 635mm	102kgs (224lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	30km (20miles) / 40km (25miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	775mm (30.5")
P13R	160kgs (350lbs)	1140mm x 635mm	105kgs (231lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	30km (20miles) / 40km (25miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	775mm (30.5")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P13	160kgs (350lbs)	1140mm x 635mm	105kgs (231lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	30km (20miles) / 40km (25miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	775mm (30.5")
P13RT	160kgs (350lbs)	1140mm x 635mm	121kgs (266lbs)	450W x 2pcs	8KPH (5MPH)	12V 50Ah x 2pcs	30km (20miles) / 40km (25miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	775mm (30.5")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

POWER MOBILITY WHEELCHAIR



MINI P14

An economical model highlights seamless frame, beautiful streamlines and a harness-free battery box.



PUZZLE P15

A model which can be disassembled in several light assemblies is a perfect travel mate.



Before Being Assembled

TEMPO HPS

An economical model offers great portability and gives a clean look.



PUZZLE S P15S

A different look of P15, coming with detachable swing-away or elevating legrests rather than a fixed foot frame.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P14	135kgs (300lbs)	940mm x 610mm	55kgs (121lbs)	450W x 2pcs	5KPH (3.1MPH)	12V 22Ah x 2pcs	20km (13miles)	230mm x 90mm (9" x 3.5")	5 Amp Off Board	720mm (28.375")
HPS	100kgs (220 lbs)	910mm x 600mm	66kgs (145lbs)	450W x 2pcs	6KPH (3.8MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	715mm (28.125")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P15	135kgs (300lbs)	940mm x 610mm	38kgs (84lbs)	180W x 2pcs	6KPH (3.8MPH)	12V 22Ah x 2pcs	20km (13miles)	320mm x 60mm (12.5" x 2.375")	5 Amp Off Board	720mm (28.375")
P15S	135kgs (300lbs)	940mm x 610mm	38kgs (84lbs)	180W x 2pcs	6KPH (3.8MPH)	12V 22Ah x 2pcs	20km (12miles)	320mm x 60mm (12.5" x 2.375")	5 Amp Off Board	720mm (28.375")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

POWER MOBILITY WHEELCHAIR



FANTASY C P17C



An Aluminum-framed pediatric power wheelchair equipped with a cute captain seat ensures a lot of fun when driving.

FANTASY R P17R



An Aluminum-framed pediatric power wheelchair equipped with a rehab seating system. This model has the seat width options of 12", 14", and 16" that the seat depth is adjustable as well. When the small angle grows up, this chair does too!



FANTASY RT P17RT




An Aluminum-framed pediatric power wheelchair equipped with a tilt-in-space rehab seating system that does good for blood circulation!



SMART RT HP9RT



A pediatric power wheelchair equipped with a tilt-in-space rehab seating system caters to all your angle's needs.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P17C	100kgs (220lbs)	700mm x 620mm	80Kgs (176 lbs)	450W x 2pcs	8KPH (5MPH)	12V 36Ah x 2pcs	30km (20miles)	320mm x 60mm (12.5" x 2.375")	5 Amp Off Board	480mm (19")
P17R	100kgs (220lbs)	700mm x 620mm	80Kgs (176lbs)	450W x 2pcs	8KPH (5MPH)	12V 36Ah x 2pcs	30km (20miles)	320mm x 60mm (12.5" x 2.375")	5 Amp Off Board	480mm (19")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P17RT	100kgs (220lbs)	700mm x 620mm	89Kgs (196lbs)	450W x 2pcs	8KPH (5MPH)	12V 36Ah x 2pcs	30km (20miles)	320mm x 60mm (12.5" x 2.375")	5 Amp Off Board	480mm (19")
HP9RT	100kgs (220lbs)	950mm x 610mm	90kgs (198lbs)	450W x 2pcs	7KPH (4.4MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	610mm (24")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

POWER MOBILITY WHEELCHAIR



VITAL C P16C

A strong heavy duty power wheelchair featured with the maximum seat width of 25" and 4-pole high power motors.



VITAL R P16R

A model that is equipped with a solid rehab seating system looks strong in structure but soft in lines.



VITAL RTL P16RTL

It equipped with tilt-in-space rehab seating system, is featured with LED indication lights and linear guideways of high stiffness, precise accuracy, and great loading capacity.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P16C	225kgs (500lbs)	1150mm x 660mm	117kgs (257lbs)	700W 4-Pole Motor x 2pcs	10KPH (6.25MPH)	12V 50Ah/62Ah x 2pcs	30km (20miles) / 40km (25miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	605mm (23.75")
P16R	225kgs (500lbs)	1100mm x 660mm	124kgs (273lbs)	700W 4-Pole Motor x 2pcs	10KPH (6.25MPH)	12V 50Ah/62Ah x 2pcs	30km (20miles) / 40km (25miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	605mm (23.75")
P16RTL	225kgs (500lbs)	1100mm x 660mm	143kgs (315lbs)	700W 4-Pole Motor x 2pcs	10KPH (6.25MPH)	12V 50Ah/62Ah x 2pcs	30km (20miles) / 40km (25miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	605mm (23.75")

ERA 607 CL P18CL




A model that features an European-styled frame and a comfortable captain seat is just great for both indoor and outdoor uses.



ERA 607 RL P18RL



A model that has rehab seating system and independent rear suspensions is featured with LED lights as standard equipments.



A model with tilt-in-space rehab seating system assures you a good ride because of its independent rear suspension system.

ERA 607 RTL P18RTL



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P18CL	205kgs (450lbs)	1150mm x 680mm	92kgs (202lbs)	450W x 2pcs	10KPH (6.25MPH)	12V 50Ah/62Ah x 2pcs	30km (20miles) / 40km (25miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	600mm (23.625")
P18RL	205kgs (450lbs)	1150mm x 680mm	94kgs (207lbs)	450W x 2pcs	10KPH (6.25MPH)	12V 50Ah/62Ah x 2pcs	30km (20miles) / 40km (25miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	600mm (23.625")
P18RTL	205kgs (450lbs)	1150mm x 680mm	104kgs (229lbs)	450W x 2pcs	10KPH (6.25MPH)	12V 50Ah/62Ah x 2pcs	30km (20miles) / 40km (25miles)	360mm x 75mm (14" x 3")	5 Amp Off Board	600mm (23.625")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.

POWER MOBILITY WHEELCHAIR



CHALLENGER CL P20CL

An outdoor power wheelchair features 16" rear wheels that are perfect for a faster stroll. Its built-in LED lighting system also greatly improves the safety level while driven at night.

CHALLENGER RTL P20RTL

A tilt-in-space rehab seating system whose speed and torque promises a fantastic driving experience.



ZEUS C P23C

- A mid-wheel drive power chair offers a small turning radius
- It has the obstacle climbing ability of 50mm



ZEUS R P23R

- A mid-wheel drive power chair incorporated with rehab seating system
- It has obstacle climbing ability of 50mm

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P20CL	225kgs (500lbs)	1220mm x 730 mm	135kgs (300lbs)	700W 4-Pole Motor x 2pcs	13KPH (8.1MPH)	12V 62Ah/80Ah x 2pcs	40km (25miles) / 50km (32miles)	410mm x 90mm (16" x 3.5")	8 Amp Off Board	850mm (33.5")
P20RTL	225kgs (500lbs)	1220mm x 730 mm	145kgs (319lbs)	700W 4-Pole Motor x 2pcs	13KPH (8.1MPH)	12V 62Ah/80Ah x 2pcs	40km (25miles) / 50km (32miles)	410mm x 90mm (16" x 3.5")	8 Amp Off Board	850mm (33.5")

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P23C	135kgs (300lbs)	1100mm x 660mm	85kgs (187lbs)	250W x 2pcs	6KPH (3.8MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	550mm (21.625")
P23R	135kgs (300lbs)	1100mm x 660mm	92kgs (202lbs)	250W x 2pcs	6KPH (3.8MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	550mm (21.625")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.