



SCOOTER	Vita	Vita S	Vita T	Bien
	S12	S12S	S12T	S15
WEIGHT CAPACITY	160kgs (350lbs)	160kgs (350lbs)	160kgs (350lbs)	160kgs (350lbs)
DIMENSION L x W	1420mm x 700mm	1420mm x 700mm	1400mm x 700mm	1400mm x 700mm
WEIGHT W/BATTERY	135kgs (300lbs)	140kgs (308lbs)	133kgs (292lbs)	138kgs (303lbs)
MOTOR TYPE	700W 4-Pole Motor	700W 4-Pole Motor	700W 4-Pole Motor	700W 4-Pole Motor
CONTROLLER	S-Drive 120Amp	S-Drive 160Amp	S-Drive 120Amp	Dynamic Rhino2 120Amp
MAX SPEED	15KPH (9.4MPH)	18KPH (11.3MPH)	15KPH (9.4MPH)	13.5KPH (8.4MPH)
BATTERY	12V 50Ah/80Ah x 2pcs	12V 50Ah/80Ah x 2pcs	12V 50Ah/80Ah x 2pcs	12V 50Ah/80Ah x 2pcs
BATTERY RANGE	30km (20miles) / 50km (31miles)	30km (20miles) / 50km (31miles)	30km (20miles) / 50km (31miles)	30km (20miles) / 50km (31miles)
DRIVE WHEEL	330mm x 120mm (13" x 5") / 360mm x 120mm (14" x 5") (Optional)	330mm x 120mm (13" x 5") / 360mm x 120mm (14" x 5") (Optional)	330mm x 120mm (13" x 5") / 360mm x 120mm (14" x 5") (Optional)	330mm x 120mm (13" x 5") / 360mm x 120mm (14" x 5") (Optional)
FRONT CASTER	330mm x 100mm (13" x 4")	330mm x 100mm (13" x 4")	330mm x 100mm (13" x 4")	330mm x 100mm (13" x 4")
REAR CASTER(ANTI-TIPPER)	None	None	None	None
CHARGER TYPE	8 Amp Off Board	8 Amp Off Board	8 Amp Off Board	8 Amp Off Board
SUSPENSIONS	Full	Full	Full	Full
TURNING RADIUS	1280mm (50")	1280mm (50")	1150mm (45")	1280mm (50")



SCOOTER	Cutie	Aviator	Aviator X	Venus
	S17	S8	S8X	S9
WEIGHT CAPACITY	135kgs (300lbs)	185kgs (400lbs)	185kgs (400lbs)	135kgs (300lbs)
DIMENSION L x W	1320mm x 670mm	1570mm x 700mm	1570mm x 800mm	1270mm x 650mm
WEIGHT W/BATTERY	135kgs (300lbs)	108kgs (238lbs)	113kgs (249lbs)	105kgs (231lbs)
MOTOR TYPE	500W 4-Pole Motor	700W 4-Pole Motor	700W 4-Pole Motor	500W 4-Pole Motor
CONTROLLER	Dynamic Rhino2 120Amp / S-Drive 120Amp	S-Drive 120Amp	S-Drive 120Amp	S-Drive 120Amp
MAX SPEED	12KPH (7.5MPH)	13KPH (8.1MPH)	13KPH (8.1MPH)	10KPH (6.25MPH)
BATTERY	12V 36Ah/50Ah x 2pcs	12V 50Ah/80Ah x 2pcs	12V 50Ah/80Ah x 2pcs	12V 36Ah/50Ah x 2pcs
BATTERY RANGE	25km (16miles) / 35km (22miles)	35km (22miles) / 50km (31miles)	35km (22miles) / 50km (31miles)	30km (20miles) / 45km (28miles)
DRIVE WHEEL	280mm x 100mm (11" x 4")	330mm x 100mm (13" x 4")	330mm x 162mm (13" x 6.5")	280mm x 100mm (11" x 4")
FRONT CASTER	280mm x 90mm (11" x 3.5")	330mm x 100mm (13" x 4")	330mm x 100mm (13" x 4")	280mm x 100mm (11" x 4")
REAR CASTER(ANTI-TIPPER)	None	None	None	50mm x 23mm (2" x 0.875")
CHARGER TYPE	5 Amp / 8 Amp Off Board	5 Amp / 8 Amp Off Board	8 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Full	Full	Full	Full
TURNING RADIUS	1200mm (47")	1160mm (45.625")	1160mm (45.625")	980mm (38")



SCOOTER	Zen	Picnic	Pixi	Royale 3
	S11	S33	S34	PT7
WEIGHT CAPACITY	135 kgs (300lbs)	115kgs (250lbs)	100kgs (220lbs)	160kgs (350lbs)
DIMENSION L x W	1070mm x 580mm	1040mm x 535 mm	995mm x 490mm	1580mm x 730mm
WEIGHT W/BATTERY	58kgs (128lbs)	40kgs (88lbs)	44kgs (97lbs)	153kgs (337lbs)
MOTOR TYPE	400W	270W	270W	1300W 4-Pole Motor
CONTROLLER	S-Drive 70Amp	S-Drive 70Amp	S-Drive 45Amp	S-Drive 160Amp
MAX SPEED	8KPH (5MPH)	6KPH (3.8MPH)	7KPH (4.4MPH)	15KPH (9.4MPH)
BATTERY	12V 22Ah x 2pcs	12V 22Ah x 2pcs	12V 12Ah x 2pcs	12V 80Ah x 2pcs
BATTERY RANGE	20km (13miles)	15km (9miles)	13km (8miles)	50km (31miles)
DRIVE WHEEL	240mm x 75mm (9.5" x 3")	200mm x 50mm (8" x 2")	200mm x 75mm (8" x 3")	410mm x 100mm (16" x 4")
FRONT CASTER	240mm x 75mm (9.5" x 3")	65mm x 25mm (2.5" x 1")	200mm x 75mm (8" x 3")	410mm x 90mm (16" x 3.5")
REAR CASTER(ANTI-TIPPER)	35mm x 17mm (1.375" x 0.625")	65mm x 25mm (2.5" x 1")	35mm x 17mm (1.375" x 0.625")	None
CHARGER TYPE	2 Amp Off Board	2 Amp Off Board	2 Amp Off Board	8 Amp Off Board
SUSPENSIONS	Rear	None	None	Full
TURNING RADIUS	980mm (38")	860mm (34")	1450mm (57")	1260mm (49.5")



SCOOTER	Royale 4	Royale 4 Plus	Royale 4X	Royale 4D
	PF7	PF7+	PF7X	PF7D
WEIGHT CAPACITY	160kgs (350lbs)	160kgs (350lbs)	160kgs (350lbs)	160kgs (350lbs)
DIMENSION L x W	1580mm x 730mm	1580mm x 730mm	1580mm x 800mm	1580mm x 730mm
WEIGHT W/BATTERY	151kgs (332lbs)	152kgs (334lbs)	158kgs (348lbs)	171kgs (376lbs)
MOTOR TYPE	1300W 4-Pole Motor	1300W 4-Pole Motor	1300W 4-Pole Motor	1300W 4-Pole Motor
CONTROLLER	S-Drive 160Amp	S-Drive 160Amp	S-Drive 160Amp	S-Drive 160Amp
MAX SPEED	15KPH (9.4MPH)	15KPH (9.4MPH)	15KPH (9.4MPH)	15KPH (9.4MPH)
BATTERY	12V 80Ah x 2pcs	12V 80Ah x 2pcs	12V 80Ah x 2pcs	12V 80Ah x 2pcs
BATTERY RANGE	50km (31miles)	50km (31miles)	50km (31miles)	50km (31miles)
DRIVE WHEEL	410mm x 100mm (16" x 4")	410mm x 100mm (16" x 4")	410mm x 150mm (16" x 6")	410mm x 100mm (16" x 4")
FRONT CASTER	360mm x 90mm (14" x 3.5")	360mm x 90mm (14" x 3.5")	360mm x 90mm (14" x 3.5")	360mm x 90mm (14" x 3.5")
REAR CASTER(ANTI-TIPPER)	None	None	None	None
CHARGER TYPE	8 Amp Off Board	8 Amp Off Board	8 Amp Off Board	8 Amp Off Board
SUSPENSIONS	Full	Full	Full	Full
TURNING RADIUS	1600mm (63")	1600mm (63")	1600mm (63")	1600mm (63")



SCOOTER	Frontier	Frontier S	Frontier X	Bolero
	PF1	PF1S	PF1X	PF2
WEIGHT CAPACITY	185kgs (400lbs)	185kgs (400lbs)	185kgs (400lbs)	160kgs (350lbs)
DIMENSION L x W	1420mm x 635mm	1420mm x 635mm	1420mm x 635mm	1310mm x 610mm
WEIGHT W/BATTERY	100kgs (220lbs)	100kgs (220lbs)	100kgs (220lbs)	90kgs (198lbs)
MOTOR TYPE	550W	550W	700W 4-Pole Motor	450W
CONTROLLER	S-Drive 120Amp	S-Drive 120Amp	S-Drive 160Amp	S-Drive 120Amp
MAX SPEED	12KPH (7.5MPH)	12KPH (7.5MPH)	13KPH (8.1MPH)	10KPH (6.25MPH)
BATTERY	12V 50Ah x 2pcs	12V 50Ah x 2pcs	12V 50Ah x 2pcs	12V 36Ah x 2pcs
BATTERY RANGE	40km (25miles)	40km (25miles)	40km (25miles)	30km (20miles)
DRIVE WHEEL	305mm x 100mm (12" x 4")	305mm x 100mm (12" x 4")	305mm x 100mm (12" x 4")	260mm x 85mm (10" x 3.375")
FRONT CASTER	305mm x 100mm (12" x 4")	305mm x 100mm (12" x 4")	305mm x 100mm (12" x 4")	260mm x 85mm (10" x 3.375")
REAR CASTER(ANTI-TIPPER)	65mm x 25mm (2.5" x 1")	65mm x 25mm (2.5" x 1")	65mm x 25mm (2.5" x 1")	50mm x 23mm (2" x 0.875")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	None	Full	Full	None
TURNING RADIUS	1160mm (45.625")	1160mm (45.625")	1160mm (45.625")	1100mm (43.25")



SCOOTER	Bolero S	Nomad 4	Nomad 3	Mystere
	PF2S	PF3	PT3	PF5
WEIGHT CAPACITY	160kgs (350lbs)	135kgs (300lbs)	135kgs (300lbs)	185kgs (400lbs)
DIMENSION L x W	1310mm x 610mm	1155mm x 600mm	1100mm x 600mm	1420mm x 630mm
WEIGHT W/BATTERY	90kgs (198lbs)	62kgs (136lbs)	54kgs (119lbs)	108kgs (238lbs)
MOTOR TYPE	450W	450W	450W	550W
CONTROLLER	S-Drive 120Amp	S-Drive 70Amp	S-Drive 70Amp	S-Drive 120Amp
MAX SPEED	7KPH (4.4MPH)	7KPH (4.4MPH)	7KPH (4.4MPH)	12KPH (7.5MPH)
BATTERY	12V 36Ah x 2pcs	12V 36Ah x 2pcs	12V 36Ah x 2pcs	12V 50Ah x 2pcs
BATTERY RANGE	30km (20miles)	30km (20miles)	30km (20miles)	40km (25miles)
DRIVE WHEEL	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	305mm x 100mm (12" x 4")
FRONT CASTER	260mm x 85mm (10" x 3.375")	240mm x 75mm (9.5" x 3")	240mm x 75mm (9.5" x 3")	305mm x 100mm (12" x 4")
REAR CASTER(ANTI-TIPPER)	50mm x 23mm (2" x 0.875")	35mm x 17mm (1.375" x 0.625")	35mm x 17mm (1.375" x 0.625")	70mm x 25mm (2.75" x 1")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Full	None	None	Full
TURNING RADIUS	1100mm (43.25")	970mm (38")	940mm (37")	1160mm (45.625")



SCOOTER	Mystere X	Mirage	Mirage C	Mirage C*
	PF5X	PF6	PF6C	PF6C*
WEIGHT CAPACITY	185kgs (400lbs)	135kgs (300lbs)	135kgs (300lbs)	135kgs (300lbs)
DIMENSION L x W	1420mm x 630mm	1220mm x 610mm	1220mm x 610mm	1220mm x 610mm
WEIGHT W/BATTERY	108kgs (238lbs)	79kgs (174lbs)	85kgs (187lbs)	86kgs (189lbs)
MOTOR TYPE	700W 4-Pole Motor	450W	450W	450W
CONTROLLER	S-Drive 160Amp	S-Drive 70Amp	S-Drive 70Amp	S-Drive 120Amp
MAX SPEED	13KPH (8.1MPH)	10KPH (6.25MPH)	10KPH (6.25MPH)	11KPH (6.8MPH)
BATTERY	12V 50Ah x 2pcs	12V 36Ah x 2pcs	12V 36Ah x 2pcs	12V 36Ah x 2pcs
BATTERY RANGE	40km (25miles)	30km (20miles)	30km (20miles)	30km (20miles)
DRIVE WHEEL	305mm x 100mm (12" x 4")	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	280mm x 100mm (11" x 4")
FRONT CASTER	305mm x 100mm (12" x 4")	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	280mm x 100mm (11" x 4")
REAR CASTER(ANTI-TIPPER)	70mm x 25mm (2.75" x 1")	50mm x 25mm (2" x 1")	50mm x 25mm (2" x 1")	50mm x 25mm (2" x 1")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Full	None	None	None
TURNING RADIUS	1160mm (45.625")	1060mm (41.75")	1060mm (41.75")	1060mm (41.75")



SCOOTER	Mirage K
	PF6K
WEIGHT CAPACITY	135kgs (300lbs)
DIMENSION L x W	1220mm x 610mm
WEIGHT W/BATTERY	87kgs (191lbs)
MOTOR TYPE	450W
CONTROLLER	S-Drive 120Amp
MAX SPEED	12KPH (7.5MPH)
BATTERY	12V 36Ah x 2pcs
BATTERY RANGE	20km (13miles)
DRIVE WHEEL	305mm x 100mm (12" x 4")
FRONT CASTER	305mm x 100mm (12" x 4")
REAR CASTER(ANTI-TIPPER)	50mm x 25mm (2" x 1")
CHARGER TYPE	5 Amp Off Board
SUSPENSIONS	None
TURNING RADIUS	1060mm (41.75")





WHEELCHAIR	CEO	CEO CL	CEO RL	CEO RTL
	P25	P25CL	P25RL	P25RTL
WEIGHT CAPACITY	160kgs (350lbs)	160kgs (350lbs)	160kgs (350lbs)	160kgs (350lbs)
DIMENSION L x W	830mm x 650mm	830mm x 650mm	830mm x 650mm	830mm x 650mm
WEIGHT W/BATTERY	135kgs (300lbs)	97kgs (213lbs)	106kgs (233lbs)	125kgs (275lbs)
MOTOR TYPE	450W 4-Pole Motor x 2pcs	450W 4-Pole Motor x 2pcs	450W 4-Pole Motor x 2pcs	450W 4-Pole Motor x 2pcs
CONTROLLER	R-net 90Amp	VR2 70Amp	VR2 70Amp	VR2 70Amp
MAX SPEED	10KPH (6.25MPH)	10KPH (6.25MPH)	10KPH (6.25MPH)	10KPH (6.25MPH)
BATTERY	12V 80Ah x 2pcs	12V 80Ah x 2pcs	12V 80Ah x 2pcs	12V 80Ah x 2pcs
BATTERY RANGE	50km (31miles)	50km (31miles)	50km (31miles)	50km (31miles)
DRIVE WHEEL	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")
FRONT CASTER	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")
REAR CASTER(ANTI-TIPPER)	75mm x 25mm (3" x 1")	75mm x 25mm (3" x 1")	75mm x 25mm (3" x 1")	75mm x 25mm (3" x 1")
CHARGER TYPE	8 Amp Off Board	8 Amp Off Board	8 Amp Off Board	8 Amp Off Board
SUSPENSIONS	Full	Full	Full	Full
TURNING RADIUS	840mm (33")	840mm (33")	840mm (33")	840mm (33")



WHEELCHAIR	Flyer	Flyer C	Flyer R	Flyer RT
	P25J	P25JC	P25JR	P25JRT
WEIGHT CAPACITY	80kgs (176lbs)	80kgs (176lbs)	80kgs (176lbs)	80kgs (176lbs)
DIMENSION L x W	920mm x 620mm	920mm x 620mm	920mm x 620mm	920mm x 620mm
WEIGHT W/BATTERY	75kgs (165lbs)	64kgs (141lbs)	65kgs (143lbs)	70kgs (154lbs)
MOTOR TYPE	250W x 2pcs	250W x 2pcs	250W x 2pcs	250W x 2pcs
CONTROLLER	R-net 90Amp	VR2 70Amp	VR2 70Amp	VR2 70Amp
MAX SPEED	9KPH (5.6MPH)	9KPH (5.6MPH)	9KPH (5.6MPH)	9KPH (5.6MPH)
BATTERY	12V 36Ah/50Ah x 2pcs	12V 36Ah/50Ah x 2pcs	12V 36Ah/50Ah x 2pcs	12V 36Ah/50Ah x 2pcs
BATTERY RANGE	25km (16miles) / 35km (22miles)	25km (16miles) / 35km (22miles)	25km (16miles) / 35km (22miles)	25km (16miles) / 35km (22miles)
DRIVE WHEEL	330mm x 75mm (13" x 3")	330mm x 75mm (13" x 3")	330mm x 75mm (13" x 3")	330mm x 75mm (13" x 3")
FRONT CASTER	230mm x 70mm (9" x 2.75")	230mm x 70mm (9" x 2.75")	230mm x 70mm (9" x 2.75")	230mm x 70mm (9" x 2.75")
REAR CASTER(ANTI-TIPPER)	60mm x 20mm (2.375" x 0.75")	60mm x 20mm (2.375" x 0.75")	60mm x 20mm (2.375" x 0.75")	60mm x 20mm (2.375" x 0.75")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Full	Full	Full	Full
TURNING RADIUS	750mm (29.5")	750mm (29.5")	750mm (29.5")	750mm (29.5")



WHEELCHAIR	KNIGHT CL	KNIGHT RL	KNIGHT RTL	Maxx
	P26CL	P26RL	P26RTL	P3DXC
WEIGHT CAPACITY	185kgs (400lbs)	185kgs (400lbs)	185kgs (400lbs)	160kgs (350lbs)
DIMENSION L x W	930mm x 600mm	930mm x 600mm	930mm x 600mm	930mm x 600mm
WEIGHT W/BATTERY	120kgs (264lbs)	145kgs (319lbs)	165kgs (363lbs)	98kgs (216lbs)
MOTOR TYPE	700W 4-Pole Motor x 2pcs	700W 4-Pole Motor x 2pcs	700W 4-Pole Motor x 2pcs	500W 4-Pole Motor x 2pcs
CONTROLLER	VR2 90Amp	VR2 90Amp	VR2 90Amp	VR2 90Amp
MAX SPEED	15KPH (9.4MPH)	15KPH (9.4MPH)	15KPH (9.4MPH)	10KPH (6.25MPH)
BATTERY	12V 80Ah x 2pcs	12V 80Ah x 2pcs	12V 80Ah x 2pcs	12V 50Ah x 2pcs
BATTERY RANGE	50km (31miles)	50km (31miles)	50km (31miles)	30km (20miles)
DRIVE WHEEL	410mm x 90mm (16" x 3.5")	410mm x 90mm (16" x 3.5")	410mm x 90mm (16" x 3.5")	360mm x 75mm (14" x 3")
FRONT CASTER	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	155mm x 50mm (6" x 2")
REAR CASTER(ANTI-TIPPER)	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")	155mm x 50mm (6" x 2")
CHARGER TYPE	8 Amp Off Board	8 Amp Off Board	8 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Drive Wheel / Front Caster	Drive Wheel / Front Caster	Drive Wheel / Front Caster	Drive Wheel / Front Caster
TURNING RADIUS	850mm (33.375")	850mm (33.375")	850mm (33.375")	500mm (19.625")



WHEELCHAIR	Maxx	Maxx	Spring	Spring C
	P3DXR	P3DXRT	HW1	HW1C
WEIGHT CAPACITY	160kgs (350lbs)	160kgs (350lbs)	135kgs (300lbs)	135kgs (300lbs)
DIMENSION L x W	930mm x 600mm	930mm x 600mm	1150mm x 700mm	1150mm x 700mm
WEIGHT W/BATTERY	107kgs (235lbs)	126kgs (277lbs)	36.6kgs (80lbs)	36.6kgs (80lbs)
MOTOR TYPE	500W 4-Pole Motor x 2pcs	500W 4-Pole Motor x 2pcs	None	None
CONTROLLER	VR2 90Amp	VR2 90Amp	None	None
MAX SPEED	10KPH (6.25MPH)	10KPH (6.25MPH)	None	None
BATTERY	12V 50Ah x 2pcs	12V 50Ah x 2pcs	None	None
BATTERY RANGE	30km (20miles)	30km (20miles)	None	None
DRIVE WHEEL	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")	610mm (24")	305mm (12")
FRONT CASTER	155mm x 50mm (6" x 2")	155mm x 50mm (6" x 2")	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")
REAR CASTER(ANTI-TIPPER)	155mm x 50mm (6" x 2")	155mm x 50mm (6" x 2")	None	None
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	None	None
SUSPENSIONS	Drive Wheel / Front Caster	Drive Wheel / Front Caster	None	None
TURNING RADIUS	500mm (19.625")	500mm (19.625")	-	-



WHEELCHAIR	Spring P	Spring J	Spring JC	Escape
	HW1P	HW1J	HW1JC	HP1
WEIGHT CAPACITY	135kgs (300lbs)	100kgs (220lbs)	100kgs (220lbs)	135kgs (300lbs)
DIMENSION L x W	1150mm x 700mm	890mm x 585mm	890mm x 585mm	1090mm x 630mm
WEIGHT W/BATTERY	58kgs (128lbs)	26kgs (57lbs)	26kgs (57lbs)	62kgs (136lbs)
MOTOR TYPE	None	None	None	200W x 2pcs
CONTROLLER	None	None	None	VSI 50Amp
MAX SPEED	None	None	None	6KPH (3.8MPH)
BATTERY	12V 1.3Ah x 2pcs	None	None	12V 36Ah x 2pcs
BATTERY RANGE	None	None	None	30km (20miles)
DRIVE WHEEL	610mm (24")	510mm (20")	305mm (12")	330mm x 50mm (13" x 2")
FRONT CASTER	200mm x 50mm (8" x 2")	178mm x 45mm (7" x 1.75")	178mm x 45mm (7" x 1.75")	200mm x 50mm (8" x 2")
REAR CASTER(ANTI-TIPPER)	None	None	None	35mm x 17mm (1.375" x 0.625")
CHARGER TYPE	None	None	None	5 Amp Off Board
SUSPENSIONS	None	None	None	None
TURNING RADIUS	-	-	-	865mm (34")



WHEELCHAIR	Escape DX	Escape LX	Escape S	Escape SX
	HP5	HP8	P12S	P12SX
WEIGHT CAPACITY	135kgs (300lbs)	135kgs (300lbs)	135kgs (300lbs)	160kgs (350lbs)
DIMENSION L x W	1100mm x 685mm	1085mm x 640mm	1120mm x 610mm	1130mm x 660mm
WEIGHT W/BATTERY	62kgs (136lbs)	66kgs (145lbs)	73kgs (161lbs)	83kgs (183lbs)
MOTOR TYPE	200W x 2pcs	200W x 2pcs	200W x 2pcs	600W x 2pcs
CONTROLLER	VSI 50Amp	VSI 90Amp	VSI 90Amp	VSI 90Amp
MAX SPEED	6KPH (3.8MPH)	6KPH (3.8MPH)	6KPH (3.8MPH)	12KPH (7.5MPH)
BATTERY	12V 36Ah x 2pcs	12V 36Ah x 2pcs	12V 36Ah/50Ah x 2pcs	12V 36Ah/50Ah x 2pcs
BATTERY RANGE	30km (20miles)	30km (20miles)	30km (20miles) / 45km (28miles)	30km (20miles) / 45km (28miles)
DRIVE WHEEL	320mm x 60mm (12.5" x 2.375")	320mm x 60mm (12.5" x 2.375")	320mm x 60mm (12.5" x 2.375")	360mm x 75mm (14" x 3")
FRONT CASTER	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")	240mm x 75mm (9.5" x 3")
REAR CASTER(ANTI-TIPPER)	35mm x 17mm (1.375" x 0.625")	35mm x 17mm (1.375" x 0.625")	35mm x 17mm (1.375" x 0.625")	35mm x 17mm (1.375" x 0.625")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	None	None	Drive Wheel	Drive Wheel
TURNING RADIUS	865mm (34")	865mm (34")	900mm (35.5")	950mm (37")



WHEELCHAIR	Escape AX	Escape HD	Escape LXHD	Maxx
	P21	HP1HD	HP8HD	P3D
WEIGHT CAPACITY	135kgs (300lbs)	180kgs (400lbs)	225kgs (500lbs)	160kgs (350lbs)
DIMENSION L x W	1030mm x 630mm	1150mm x 750mm	1170mm x 910mm	1080mm x 680mm
WEIGHT W/BATTERY	80kgs (176lbs)	82kgs (180lbs)	85kgs (187lbs)	107kgs (235lbs)
MOTOR TYPE	450W x 2pcs	700W 4-Pole Motor x 2pcs	700W 4-Pole Motor x 2pcs	450W x 2pcs
CONTROLLER	VSI 50Amp	VR2 90Amp	VR2 90Amp	VR2 50Amp
MAX SPEED	6KPH (3.8MPH)	10KPH (6.25MPH)	10KPH (6.25MPH)	8KPH (5MPH)
BATTERY	12V 50Ah x 2pcs	12V 50Ah x 2pcs	12V 62Ah x 2pcs	12V 50Ah x 2pcs
BATTERY RANGE	45km (28miles)	20km (13miles)	30km (20miles)	30km (20miles)
DRIVE WHEEL	320mm x 60mm (12.5" x 2.375")	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")	305mm x 100mm (12" x 4")
FRONT CASTER	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")	155mm x 50mm (6" x 2")
REAR CASTER(ANTI-TIPPER)	35mm x 17mm (1.375" x 0.625")	35mm x 17mm (1.375" x 0.625")	35mm x 17mm (1.375" x 0.625")	155mm x 50mm (6" x 2")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	8 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Drive Wheel	None	None	Drive Wheel / Front Caster
TURNING RADIUS	850mm (33.5")	1000mm (39")	1140mm (45")	550mm (21.625")



WHEELCHAIR	Maxx R	Maxx RT	Rumba	Rumba S
	P3DR	P3DRT	HP3	HP4
WEIGHT CAPACITY	160kgs (350lbs)	160kgs (350lbs)	135kgs (300lbs)	135kgs (300lbs)
DIMENSION L x W	1080mm x 680mm	1080mm x 680mm	1000mm x 610mm	1100mm x 630mm
WEIGHT W/BATTERY	107kgs (235lbs)	107kgs (235lbs)	98kgs (216lbs)	83kgs (183lbs)
MOTOR TYPE	450W x 2pcs	450W x 2pcs	450W x 2pcs	450W x 2pcs
CONTROLLER	VR2 50Amp	VR2 50Amp	VSI 50Amp	VSI 50Amp
MAX SPEED	8KPH (5MPH)	8KPH (5MPH)	8KPH (5MPH)	7KPH (4.4MPH)
BATTERY	12V 50Ah x 2pcs	12V 50Ah x 2pcs	12V 36Ah/50Ah x 2pcs	12V 36Ah/50Ah x 2pcs
BATTERY RANGE	30km (20miles)	30km (20miles)	30km (20miles) / 45km (28miles)	30km (20miles) / 45km (28miles)
DRIVE WHEEL	330mm x 100mm (13" x 4")	330mm x 100mm (13" x 4")	330mm x 100mm (13" x 4")	260mm x 85mm (10" x 3.375")
FRONT CASTER	155mm x 50mm (6" x 2")	155mm x 50mm (6" x 2")	240mm x 75mm (9.5" x 3")	200mm x 50mm (8" x 2")
REAR CASTER(ANTI-TIPPER)	155mm x 50mm (6" x 2")	155mm x 50mm (6" x 2")	50mm x 23mm (2" x 0.875")	50mm x 23mm (2" x 0.875")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Drive Wheel / Front Caster	Drive Wheel / Front Caster	None	None
TURNING RADIUS	550mm (21.625")	550mm (21.625")	825mm (32.5")	785mm (31")



WHEELCHAIR	Rumba SF	Rumba SR	Tiara	Tiara R
	P4F	P4R	P4AS	P4ASR
WEIGHT CAPACITY	135kgs (300lbs)	135kgs (300lbs)	135kgs (300lbs)	135kgs (300lbs)
DIMENSION L x W	1110mm x 610mm	1040mm x 610mm	1100mm x 660mm	1100mm x 660mm
WEIGHT W/BATTERY	93kgs (205lbs)	93kgs (205lbs)	90kgs (198lbs)	94kgs (207lbs)
MOTOR TYPE	450W x 2pcs	450W x 2pcs	450W x 2pcs	450W x 2pcs
CONTROLLER	VSI 50Amp	VSI 50Amp	VR2 60Amp	VR2 60Amp
MAX SPEED	7KPH (4.4MPH)	7KPH (4.4MPH)	7KPH (4.4MPH)	7KPH (4.4MPH)
BATTERY	12V 36Ah/50Ah x 2pcs	12V 36Ah x 2pcs	12V 36Ah x 2pcs	12V 36Ah x 2pcs
BATTERY RANGE	30km (20miles) / 45km (28miles)	30km (20miles) / 45km (28miles)	30km (20miles)	30km (20miles)
DRIVE WHEEL	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")
FRONT CASTER	100mm x 35mm (4" x 1.375")	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")
REAR CASTER(ANTI-TIPPER)	200mm x 50mm (8" x 2")	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	None	None	ANTI-TIPPER WHEEL	ANTI-TIPPER WHEEL
TURNING RADIUS	530mm (20.875")	755mm (29.5")	790mm (31")	860mm (34")



WHEELCHAIR	Tango	Tango S	Allure	Allure S
	TC1	TC1S	HP6	HP6S
WEIGHT CAPACITY	135kgs (300lbs)	135kgs (300lbs)	160kgs (350lbs)	160kgs (350lbs)
DIMENSION L x W	1160mm x 590mm	1160mm x 590mm	1110mm x 620mm	1110mm x 620mm
WEIGHT W/BATTERY	77kgs (169lbs)	77kgs (169lbs)	97kgs (213lbs)	97kgs (213lbs)
MOTOR TYPE	450W x 2pcs	450W x 2pcs	450W x 2pcs	450W x 2pcs
CONTROLLER	S-Drive 70Amp	VSI 50Amp	VSI 50Amp	VSI 50Amp
MAX SPEED	6KPH (3.8MPH)	6KPH (3.8MPH)	8KPH (5MPH)	8KPH (5MPH)
BATTERY	12V 36Ah x 2pcs	12V 36Ah x 2pcs	12V 50Ah x 2pcs	12V 50Ah x 2pcs
BATTERY RANGE	30km (20miles)	30km (20miles)	40km (25miles)	40km (25miles)
DRIVE WHEEL	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")	330mm x 100mm (13" x 4")	330mm x 100mm (13" x 4")
FRONT CASTER	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")
REAR CASTER(ANTI-TIPPER)	65mm x 25mm (2.5" x 1")	65mm x 25mm (2.5" x 1")	155mm x 35mm (6" x 1.375")	155mm x 35mm (6" x 1.375")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	None	None	Rear Caster	Front Caster and Rear Caster
TURNING RADIUS	750mm (29.5")	750mm (29.5")	660mm (26")	660mm (26")



WHEELCHAIR	Allure R	Allure RT	Sahara R	Sahara RT
	HP6R	HP6RT	HP7R	HP7RT
WEIGHT CAPACITY	160kgs (350lbs)	160kgs (350lbs)	205kgs (450lbs)	205kgs (450lbs)
DIMENSION L x W	1170mm x 620mm	1170mm x 620mm	1170mm x 635mm	1170mm x 635mm
WEIGHT W/BATTERY	101kgs (222lbs)	114kgs (251lbs)	113kgs (249lbs)	135kgs (300lbs)
MOTOR TYPE	450W x 2pcs	450W x 2pcs	700W 4-Pole Motor x 2pcs	700W 4-Pole Motor x 2pcs
CONTROLLER	VSI 50Amp	VSI 50Amp+A1	VR2 90Amp	VR2 90Amp+A1
MAX SPEED	8KPH (5MPH)	8KPH (5MPH)	10KPH (6.25MPH)	10KPH (6.25MPH)
BATTERY	12V 50Ah x 2pcs	12V 50Ah/62Ah x 2pcs	12V 50Ah/62Ah x 2pcs	12V 50Ah/62Ah x 2pcs
BATTERY RANGE	40km (25miles)	40km (25miles)	20km (13miles) / 30km (20miles)	20km (13miles) / 30km (20miles)
DRIVE WHEEL	330mm x 100mm (13" x 4")	330mm x 100mm (13" x 4")	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")
FRONT CASTER	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")	240mm x 75mm (9.5" x 3")	240mm x 75mm (9.5" x 3")
REAR CASTER(ANTI-TIPPER)	155mm x 35mm (6" x 1.375")	155mm x 35mm (6" x 1.375")	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Rear Caster	Rear Caster	Full	Full
TURNING RADIUS	660mm (26")	660mm (26")	865mm (34")	865mm (34")



WHEELCHAIR	Sahara K	Sahara KX	Vision	Vision C
	HP7K	HP7KX	P13	P13C
WEIGHT CAPACITY	205kgs (450lbs)	205kgs (450lbs)	160kgs (350lbs)	160kgs (350lbs)
DIMENSION L x W	1170mm x 830mm	1170mm x 830mm	1140mm x 635mm	1120mm x 635mm
WEIGHT W/BATTERY	128kgs (282lbs)	128kgs (282lbs)	105kgs (231lbs)	102kgs (224lbs)
MOTOR TYPE	700W 4-Pole Motor x 2pcs	700W 4-Pole Motor x 2pcs	450W x 2pcs	450W x 2pcs
CONTROLLER	VR2 90Amp	VR2 90Amp	VR2 90Amp+A1	VR2 60Amp
MAX SPEED	10KPH (6.25MPH)	10KPH (6.25MPH)	8KPH (5MPH)	8KPH (5MPH)
BATTERY	12V 50Ah/62Ah x 2pcs	12V 50Ah/62Ah x 2pcs	12V 50Ah x 2pcs	12V 50Ah x 2pcs
BATTERY RANGE	20km (13miles) / 30km (20miles)	30km (20miles)	30km (20miles) / 40km (25miles)	30km (20miles) / 40km (25miles)
DRIVE WHEEL	360mm x 160mm (14" x 6.5")	360mm x 160mm (14" x 6.5")	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")
FRONT CASTER	240mm x 75mm (9.5" x 3")	240mm x 75mm (9.5" x 3")	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")
REAR CASTER(ANTI-TIPPER)	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")
CHARGER TYPE	8 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Full	Full	Anti-tipper Wheel	Anti-tipper Wheel
TURNING RADIUS	865mm (34")	865mm (34")	775mm (30.5")	775mm (30.5")



WHEELCHAIR	Vision R	Vision RT	Mini	Tempo
	P13R	P13RT	P14	HPS
WEIGHT CAPACITY	160kgs (350lbs)	160kgs (350lbs)	135kgs (300lbs)	100kgs (220 lbs)
DIMENSION L x W	1140mm x 635mm	1140mm x 635mm	940mm x 610mm	910mm x 600mm
WEIGHT W/BATTERY	105kgs (231lbs)	121kgs (266lbs)	55kgs (121lbs)	66kgs (145lbs)
MOTOR TYPE	450W x 2pcs	450W x 2pcs	450W x 2pcs	450W x 2pcs
CONTROLLER	VR2 60Amp	VR2 60Amp+A1	VR2 50Amp	VSI 50Amp
MAX SPEED	8KPH (5MPH)	8KPH (5MPH)	5KPH (3.1MPH)	6KPH (3.8MPH)
BATTERY	12V 50Ah x 2pcs	12V 22Ah x 2pcs	12V 22Ah x 2pcs	12V 36Ah x 2pcs
BATTERY RANGE	30km (20miles) / 40km (25miles)	30km (20miles) / 40km (25miles)	20km (13miles)	30km (20miles)
DRIVE WHEEL	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")	230mm x 90mm (9" x 3.5")	260mm x 85mm (10" x 3.375")
FRONT CASTER	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")	180mm x 45mm (7" x 1.75")	200mm x 50mm (8" x 2")
REAR CASTER(ANTI-TIPPER)	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")	65mm x 25mm (2.5" x 1")	70mm x 25mm (2.75" x 1")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Anti-tipper Wheel	Anti-tipper Wheel	None	None
TURNING RADIUS	775mm (30.5")	775mm (30.5")	720mm (28.375")	715mm (28.125")



WHEELCHAIR	Puzzle	Puzzle S	Vital C	Vital R
	P15	P15S	P16C	P16R
WEIGHT CAPACITY	135kgs (300lbs)	135kgs (300lbs)	225kgs (500lbs)	225kgs (500lbs)
DIMENSION L x W	940mm x 610mm	940mm x 610mm	1150mm x 660mm	1100mm x 660mm
WEIGHT W/BATTERY	38kgs (84lbs)	38kgs (84lbs)	117kgs (257lbs)	124kgs (273lbs)
MOTOR TYPE	180W x 2pcs	180W x 2pcs	700W 4-Pole Motor x 2pcs	700W 4-Pole Motor x 2pcs
CONTROLLER	VR2 50Amp	VR2 50Amp	VR2 90Amp	VR2 90Amp
MAX SPEED	6KPH (3.8MPH)	6KPH (3.8MPH)	10KPH (6.25MPH)	10KPH (6.25MPH)
BATTERY	12V 22Ah x 2pcs	12V 22Ah x 2pcs	12V 50Ah/62Ah x 2pcs	12V 50Ah/62Ah x 2pcs
BATTERY RANGE	20km (13miles)	20km (12miles)	30km (20miles) / 40km (25miles)	30km (20miles) / 40km (25miles)
DRIVE WHEEL	320mm x 60mm (12.5" x 2.375")	320mm x 60mm (12.5" x 2.375")	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")
FRONT CASTER	180mm x 45mm (7" x 1.75")	180mm x 45mm (7" x 1.75")	240mm x 75mm (9.5" x 3")	240mm x 75mm (9.5" x 3")
REAR CASTER(ANTI-TIPPER)	35mm x 17mm (1.375" x 0.625")	35mm x 17mm (1.375" x 0.625")	130mm x 25mm (5" x 1")	130mm x 25mm (5" x 1")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	None	None	Drive Wheel / Rear Caster	Drive Wheel / Rear Caster
TURNING RADIUS	720mm (28.375")	720mm (28.375")	605mm (23.75")	605mm (23.75")



WHEELCHAIR	Vital RTL	Fantasy C	Fantasy R	Fantasy RT
	P16RTL	P17C	P17R	P17RT
WEIGHT CAPACITY	225kgs (500lbs)	100kgs (220lbs)	100kgs (220lbs)	100kgs (220lbs)
DIMENSION L x W	1100mm x 660mm	700mm x 620mm	700mm x 620mm	700mm x 620mm
WEIGHT W/BATTERY	143kgs (315lbs)	80Kgs (176 lbs)	80Kgs (176lbs)	89Kgs (196lbs)
MOTOR TYPE	700W 4-Pole Motor x 2pcs	450W x 2pcs	450W x 2pcs	450W x 2pcs
CONTROLLER	VR2 90Amp	VR2 60Amp	VR2 60Amp	VR2 60Amp
MAX SPEED	10KPH (6.25MPH)	8KPH (5MPH)	8KPH (5MPH)	8KPH (5MPH)
BATTERY	12V 50Ah/62Ah x 2pcs	12V 36Ah x 2pcs	12V 36Ah x 2pcs	12V 36Ah x 2pcs
BATTERY RANGE	30km (20miles) / 40km (25miles)	30km (20miles)	30km (20miles)	30km (20miles)
DRIVE WHEEL	360mm x 75mm (14" x 3")	320mm x 60mm (12.5" x 2.375")	320mm x 60mm (12.5" x 2.375")	320mm x 60mm (12.5" x 2.375")
FRONT CASTER	240mm x 75mm (9.5" x 3")	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")	200mm x 50mm (8" x 2")
REAR CASTER(ANTI-TIPPER)	130mm x 25mm (5" x 1")	70mm x 25mm (2.75" x 1")	70mm x 25mm (2.75" x 1")	70mm x 25mm (2.75" x 1")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Drive Wheel / Rear Caster	Drive Wheel	Drive Wheel	Drive Wheel
TURNING RADIUS	605mm (23.75")	480mm (19")	480mm (19")	480mm (19")




WHEELCHAIR	Smart RT	Era 607 CL	Era 607 RL	Era 607 RTL
	HP9RT	P18CL	P18RL	P18RTL
WEIGHT CAPACITY	100kgs (220lbs)	205kgs (450lbs)	205kgs (450lbs)	205kgs (450lbs)
DIMENSION L x W	950mm x 610mm	1150mm x 680mm	1150mm x 680mm	1150mm x 680mm
WEIGHT W/BATTERY	90kgs (198lbs)	92kgs (202lbs)	94kgs (207lbs)	104kgs (229lbs)
MOTOR TYPE	450W x 2pcs	450W x 2pcs	450W x 2pcs	450W x 2pcs
CONTROLLER	VSI 50Amp	VR2 70Amp	VR2 70Amp	VR2 70Amp+A1
MAX SPEED	7KPH (4.4MPH)	10KPH (6.25MPH)	10KPH (6.25MPH)	10KPH (6.25MPH)
BATTERY	12V 36Ah x 2pcs	12V 50Ah/62Ah x 2pcs	12V 50Ah/62Ah x 2pcs	12V 50Ah/62Ah x 2pcs
BATTERY RANGE	30km (20miles)	30km (20miles) / 40km (25miles)	30km (20miles) / 40km (25miles)	30km (20miles) / 40km (25miles)
DRIVE WHEEL	260mm x 85mm (10" x 3.375")	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")	360mm x 75mm (14" x 3")
FRONT CASTER	100mm x 35mm (4" x 1.375")	240mm x 75mm (9.5" x 3")	240mm x 75mm (9.5" x 3")	240mm x 75mm (9.5" x 3")
REAR CASTER(ANTI-TIPPER)	200mm x 50mm (8" x 2")	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")
CHARGER TYPE	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Anti-tipper Wheel	Drive Wheel	Drive Wheel	Drive Wheel
TURNING RADIUS	610mm (24")	600mm (23.625")	600mm (23.625")	600mm (23.625")



WHEELCHAIR	Challenger CL	Challenger RTL	Zeus C	Zeus R
	P20CL	P20RTL	P23C	P23R
WEIGHT CAPACITY	225kgs (500lbs)	225kgs (500lbs)	135kgs (300lbs)	135kgs (300lbs)
DIMENSION L x W	1220mm x 730 mm	1220mm x 730 mm	1100mm x 660mm	1100mm x 660mm
WEIGHT W/BATTERY	135kgs (300lbs)	145kgs (319lbs)	85kgs (187lbs)	92kgs (202lbs)
MOTOR TYPE	700W 4-Pole Motor x 2pcs	700W 4-Pole Motor x 2pcs	250W x 2pcs	250W x 2pcs
CONTROLLER	VR2 90Amp	VR2 90Amp	VR2 50Amp	VR2 50Amp
MAX SPEED	13KPH (8.1MPH)	13KPH (8.1MPH)	6KPH (3.8MPH)	6KPH (3.8MPH)
BATTERY	12V 62Ah/80Ah x 2pcs	12V 62Ah/80Ah x 2pcs	12V 36Ah x 2pcs	12V 36Ah x 2pcs
BATTERY RANGE	40km (25miles) / 50km (32miles)	40km (25miles) / 50km (32miles)	30km (20miles)	30km (20miles)
DRIVE WHEEL	410mm x 90mm (16" x 3.5")	410mm x 90mm (16" x 3.5")	260mm x 85mm (10" x 3.375")	260mm x 85mm (10" x 3.375")
FRONT CASTER	240mm x 75mm (9.5" x 3")	240mm x 75mm (9.5" x 3")	155mm x 50mm (6" x 2")	155mm x 50mm (6" x 2")
REAR CASTER(ANTI-TIPPER)	100mm x 35mm (4" x 1.375")	100mm x 35mm (4" x 1.375")	155mm x 50mm (6" x 2")	155mm x 50mm (6" x 2")
CHARGER TYPE	8 Amp Off Board	8 Amp Off Board	5 Amp Off Board	5 Amp Off Board
SUSPENSIONS	Full	Full	Drive Wheel	Drive Wheel
TURNING RADIUS	850mm (33.5")	850mm (33.5")	550mm (21.625")	550mm (21.625")



R1
ONE-PIECE RAMP



R2
TWO-PIECE RAMP



Load Capacity : 350kgs(770 lbs)

1



Item No. **94209002**
Size : 280mm x 90mm

2



Item No. **94219028**
Size : 330mm x 100mm

3



Item No. **93127137**
Size : 360mm x 75mm

4



Item No. **93426001**
Size : 260mm x 80mm

5



Item No. **94219031**
Size : 330mm x 90mm

6



Item No. **74013069**
Size : 340mm x 140mm

7



Item No. **94189008**
Size : 360mm x 90mm

8



Item No. **94219024**
Size : 330mm x 115mm