

FRONTIER PF1

A big-sized scooter with soft streamlines makes a true fashion statement!



FRONTIER S PF1S

A big-sized scooter provides great comfort due to a 4-wheel independent suspension system incorporated.



Rear Suspension System

FRONTIER X PF1X

A big-sized scooter modified from PF1, coming with a 700W 4-pole motor and a 160AMP controller, is ready to surprise you with its speed and torque.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
PF1	185kgs (400lbs)	1420mm x 635mm	100kgs (220lbs)	550W	12KPH (7.5MPH)	12V 50Ah x 2pcs	40km (25miles)	305mm x 100mm (12" x 4")	5 Amp Off Board	1160mm (45.625")
PF1S	185kgs (400lbs)	1420mm x 635mm	100kgs (220lbs)	550W	12KPH (7.5MPH)	12V 50Ah x 2pcs	40km (25miles)	305mm x 100mm (12" x 4")	5 Amp Off Board	1160mm (45.625")
PF1X	185kgs (400lbs)	1420mm x 635mm	100kgs (220lbs)	700W 4-Pole Motor	13KPH (8.1MPH)	12V 50Ah x 2pcs	40km (25miles)	305mm x 100mm (12" x 4")	5 Amp Off Board	1160mm (45.625")

BOLERO PF2

A medium-sized scooter that has a great shape and is easy to drive.



BOLERO S PF2S

A medium-sized scooter with 4-wheel independent suspension system offers a super comfortable ride.



Rear Suspension System



Delta Handle Bar




Steel Bumper Design (Optional)

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
PF2	160kgs (350lbs)	1310mm x 610mm	90kgs (198lbs)	450W	10KPH (6.25MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	1100mm (43.25")
PF2S	160kgs (350lbs)	1310mm x 610mm	90kgs (198lbs)	450W	10KPH (6.25MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	1100mm (43.25")



NOMAD 4 PF3
A medium-sized 4-wheeled scooter provides great stability and comfort.



NOMAD 3 PT3
A medium-sized 3-wheeled scooter that provides great stability and comfort.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
PF3	135kgs (300lbs)	1155mm x 600mm	62kgs (136lbs)	450W	7KPH (4.4MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	970mm (38")
PT3	135kgs (300lbs)	1100mm x 600mm	54kgs (119lbs)	450W	7KPH (4.4MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	940mm (37")

MYSTERE PF5



A big-sized scooter that is great in shape and reliable in look, highlights a large digital control panel, a 4-wheel independent suspension system, Delta-handle bar, unlimited tiller ram angle adjustments and offers lots of comfort and stability.



Luxurious Appearance



LED Control Panel



Compact Rear Shape

MYSTERE X PF5X



A big-sized scooter modified from PF5, coming with a 700W 4-pole motor and a 160AMP controller, copes well with hilly places with its speed and torque.

Optional Accessories



Canopy



Oxygen Holder



Crutch Holder



Rear Basket



Rear Bag

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
PF5	185kgs (400lbs)	1420mm x 630mm	108kgs (238lbs)	550W	12KPH (7.5MPH)	12V 50Ah x 2pcs	40km (25miles)	305mm x 100mm (12" x 4")	5 Amp Off Board	1160mm (45.625")
PF5X	185kgs (400lbs)	1420mm x 630mm	108kgs (238lbs)	700W 4-Pole Motor	13KPH (8.1MPH)	12V 50Ah x 2pcs	40km (25miles)	305mm x 100mm (12" x 4")	5 Amp Off Board	1160mm (45.625")

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.



MIRAGE PF6



A medium-sized scooter that is delicately designed offers unlimited angle adjustments on account of the unique tiller ram mechanism.

MIRAGE C⁺ PF6C⁺



A medium-sized scooter equipped with a captain seat and 11" black wheel.



MIRAGE C PF6C



A medium-sized scooter modified from PF6, coming with a Captain seat instead of a plastic one.



MIRAGE K PF6K

A medium-sized scooter modified from PF6 is equipped with a Captain seat and 12" wheels which assure a faster ride.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
PF6	135kgs (300lbs)	1220mm x 610mm	79kgs (174lbs)	450W	10KPH (6.25MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	1060mm (41.75")
PF6C	135kgs (300lbs)	1220mm x 610mm	85kgs (187lbs)	450W	10KPH (6.25MPH)	12V 36Ah x 2pcs	30km (20miles)	260mm x 85mm (10" x 3.375")	5 Amp Off Board	1060mm (41.75")

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
PF6C ⁺	135kgs (300lbs)	1220mm x 610mm	86kgs (189lbs)	450W	11KPH (6.8MPH)	12V 36Ah x 2pcs	30km (20miles)	280mm x 100mm (11" x 4")	5 Amp Off Board	1060mm (41.75")
PF6K	135kgs (300lbs)	1220mm x 610mm	87kgs (191lbs)	450W	12KPH (7.5MPH)	12V 36Ah x 2pcs	20km (13miles)	305mm x 100mm (12" x 4")	5 Amp Off Board	1060mm (41.75")